

N8-acetylspermidine oxidase (propane-1,3-diamine-forming)

Cat. No. EXWM-1542

Lot. No. (See product label)

Introduction

Description Also active with N1-acetylspermine, weak activity with N1,N12-diacetylspermine.

No activity with diaminopropane, putrescine, cadaverine, diaminohexane, norspermidine, spermine and spermidine. Absence of monoamine oxidase (EC 1.4.3.4) activity. Differs in specificity from EC 1.5.3.13 (N1-acetylpolyamine oxidase), EC 1.5.3.14 (polyamine oxidase (propane-1,3-diamine-forming)), EC 1.5.3.16 (spermine oxidase) and EC 1.5.3.17 (non-specific polyamine oxidase).

Product Information

Form Liquid or lyophilized powder

EC Number EC 1.5.3.15

Reaction N8-acetylspermidine + O2 + H2O = propane-1,3-diamine + 4-acetamidobutanal +

H2O2

Notes This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C \sim -80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1