

rhodotorulapepsin

Cat. No. EXWM-4264

Lot. No. (See product label)

Introduction

Description From the imperfect yeast Rhodotorula glutinis. Somewhat similar enzymes have

been isolated from the imperfect yeast-like organism Cladosporium sp. and the

imperfect fungus Paecilomyces varioti.

Synonyms Rhodotorula aspartic proteinase; Cladosporium acid protease; Cladosporium acid

proteinase; Paecilomyces proteinase; Cladosporium aspartic proteinase; Paecilomyces proteinase; Rhodotorula glutinis aspartic proteinase; Rhodotorula glutinis acid proteinase; Rhodotorula glutinis aspartic proteinase II; Rhodotorula

acid proteinase

Product Information

Form Liquid or lyophilized powder

EC Number EC 3.4.23.26

CAS No. 37259-59-9

Reaction Specificity similar to that of pepsin A. Cleaves Z-Lys+Ala-Ala and activates

trypsinogen

Notes This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C \sim -80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1