

## **lupeol** synthase

Cat. No. EXWM-5583

Lot. No. (See product label)

## Introduction

**Description** Also forms some β-amyrin. The recombinant enzyme from Arabidopsis thaliana gives a 1:1 mixture of

lupeol and lupan-3 $\beta$ ,20-diol with small amounts of  $\beta$ -amyrin, germanicol, taraxasterol and  $\psi$ -taraxasterol.

See EC 4.2.1.128 (lupan-3β,20-diol synthase).

**Synonyms** LUPI; BPW; RcLUS

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 5.4.99.41

**Reaction** (3S)-2,3-epoxy-2,3-dihydrosqualene = lupeol

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

## Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1