

Protease Inhibitor Mix HP PLUS

Cat. No. INBZ-012

Lot. No. (See product label)

Introduction

Applications Mixture of 6 protease inhibitors with broad spectrum of activity for protease. It is especially

recommended for purification of polyHis-tagged proteins and for other applications, where metal-

chelators should be avoided.

Product Information

Formulation In 20 mM Tris-HCl buffer, pH 8.0, 10% glycerol, 100 mM NaCl.

Targets Aspartate-; cysteine-; serine proteases as well as aminopeptidases; Thermolysin and other microbial

metalloendoproteases

Storage and Shipping Information

Storage -15 °C to -25 °C

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1