

## Protease Inhibitor Mix HP PLUS

Cat. No. INBZ-012

Lot. No. (See product label)

### Introduction

#### Applications

Mixture of 6 protease inhibitors with broad spectrum of activity for protease. It is especially recommended for purification of polyHis-tagged proteins and for other applications, where metal-chelators should be avoided.

### Product Information

#### Formulation

In 20 mM Tris-HCl buffer, pH 8.0, 10% glycerol, 100 mM NaCl.

#### Targets

Aspartate-; cysteine-; serine proteases as well as aminopeptidases; Thermolysin and other microbial metalloendoproteases

### Storage and Shipping Information

#### Storage

-15 °C to -25 °C