

Oxidase Reagent

Cat. No. NATZ-040

Lot. No. (See product label)

Introduction

Description

Presence of cytochrome oxidase can be detected through the use of an oxidase reagent, which acts as an electron donor to cytochrome oxidase. If the bacteria oxidize the reagent, it turns dark purple due to the formation of indophenol blue indicating a positive reaction. This test is fit for screening species of *Neisseria*, *Aeromonas*, *Vibrio*, *Campylobacter* and *Pseudomonas*, which show a positive reaction, or for *Enterobacteriaceae*, which show a negative reaction.

Product Information

Form

liquid