

## Collagen from Achilles Tendon

Cat. No. CEFX-003

Lot. No. (See product label)

### Introduction

**Description** All non-collagenous tissues of fresh bovine achilles tendon were removed and cut into small pieces. Extract at 0°C for 6 days and replace the Na<sub>2</sub>HP0<sub>4</sub> solution several times to remove soluble proteins. Then extracted with 25% KCl solution at 0°C for 6 days to remove mucopolysaccharides. Finally, the residue was washed thoroughly with water, dehydrated with anhydrous ethanol, and air-dried. Insoluble; Nondenaturing Particle size: Shredded (Approximately 1 cm fibers)

### Product Information

<b>Species</b>	Bovine
<b>Source</b>	Achilles Tendon
<b>Form</b>	Powder