

Phosphatase Inhibitor Cocktail B

Cat. No. CEI-0786

Lot. No. (See product label)

Introduction

Description Reconstitute in 1ml water to make a solution with 200mM imidazole, 100mM

sodium fluoride, 115mM sodium molybdate, 100mM sodium orthovanadate and 400mM sodium tartrate. For the inhibition of protein tyrosine phosphatases, acid-

and alkaline-phosphatases.

Product Information

Appearance Lyophilized solid

Targets Phosphatase

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1