

Protease Inhibitor Cocktail II

Cat. No. CEI-0796

Lot. No. (See product label)

Introduction

Description Inhibits aspartic, cysteine, serine, and metalloproteases, as well as aminopeptidases. Recommended for use with bacterial cell extracts. Reconstitute each vial with 1ml DMSO and 4ml H₂O to obtain a 5ml stock solution containing 20 mM AEBSF.HCl, 1.7mM bestatin, 200uM E-64, 85mM EDTA, and 2mM pepstatin A.

Product Information

Appearance Lyophilized solid

Targets Protease