

Fast Red TR/Naphthol AS-MX Tablets

Cat. No. CSUB-0024

Lot. No. (See product label)

Introduction

Applications Developed for use as a precipitating substrate for the detection of alkaline phosphatase activity in immunochemistry. Primary uses are in immunohistology, immunoblotting, or dot blotting. Fast Red TR/Naphthol AS-MX tablets contain 0.6 mM levamisole to block endogenous alkaline phosphatase activity. Fast Red TR/Naphthol AS-MX Tablets require no additional steps or buffers to prepare an active substrate solution.

Product Information

Form tablet

Substrates Alkaline Phosphatase

Usage and Packaging

Reconstitution Each tablet set dissolved in 10 mL of deionized water yields a ready-to-use buffered solution.

Storage and Shipping Information

Storage -20°C