

ABTS Tablets

Cat. No. CSUB-0407

Lot. No. (See product label)

Introduction

Applications ABTS Tablets are the ideal substrate for enzyme immunoassays with horseradish peroxidase (HRP) as a marker enzyme.

Product Information

Solubility water: soluble

Substrates Peroxidase Substrates

Usage and Packaging

Preparation Working solution: Dissolve one ABTS tablet (5 mg) in 5 ml of the ABTS buffer solution. This solution has a light-green color. In microtiter plates, it is colorless A405 nm/1 cm should be <0.16. Work as clean as possible. Heavy metal ions or traces of peroxidase can produce a dark green ABTS solution. Such solutions are not suitable for the assay.

Storage and Shipping Information

Storage 2-8°C

Shipping Conditions wet ice