

## DAB with Metal Enhancer

Cat. No. CSUB-0435

Lot. No. (See product label)

### Introduction

**Applications** FAST DAB (3,3'-Diaminobenzidine tetrahydrochloride) with Metal Enhancer Tablets have been developed for use in immunohistology as a precipitating substrate for the localization of peroxidase activity. The DAB reaction is enhanced by the addition of CoCl<sub>2</sub> and produces a more intense stain which is resistant to alcohol.

**Synonyms** 3, 3'-Diaminobenzidine Peroxidase Substrate; DAB Peroxidase Substrate

### Product Information

**Form** tablet

**Substrates** Peroxidase Substrates

### Usage and Packaging

**Reconstitution** Each tablet set dissolved in 5 mL deionized water yields a ready-to-use buffered solution containing DAB with a metal enhancer.

### Storage and Shipping Information

**Storage** -20°C