

Pefachrome Try

Cat. No. CSUB-0595

Lot. No. (See product label)

Introduction

Applications Pefachrome substrates are used for chromogenic assays of various serine proteinases. Pefachrome Try is

used for the detection of serine proteases such as trypsin.

Synonyms Bz-Val-Gly-Arg-4-nitroanilide acetate; N-Benzoyl-L-valyl-glycyl-L-arginine-4-nitroanilide acetate

Product Information

CAS No. 161572-28-7

Molecular C26H34N8O6 ⋅ C2H4O2

Formula

Molecular 614.65

Weight

Substrates Trypsin Substrates

Storage and Shipping Information

Storage 2-8°C

Tel: 1-631-562-8517 1-516-512-3133

Email: info@creative-enzymes.com

1/1