

tellurite methyltransferase

Cat. No. EXWM-1870

Lot. No. (See product label)

Introduction

Description The enzyme is involved in the detoxification of tellurite. It can also methylate selenite and selenium dioxide.

Synonyms TehB

Product Information

Form Liquid or lyophilized powder

EC Number EC 2.1.1.265

Reaction S-adenosyl-L-methionine + tellurite = S-adenosyl-L-homocysteine + methanetelluronate

Notes This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.