

pseudaminic acid synthase

Cat. No. EXWM-2836

Lot. No. (See product label)

Introduction

Description The enzyme requires a divalent metal ion, the highest activity values are observed in the presence of

Mn2+ and Co2+ (10 mM).

Synonyms Psel; NeuB3

Product Information

Form Liquid or lyophilized powder

EC Number EC 2.5.1.97

Reaction phosphoenolpyruvate + 2,4-bis(acetylamino)-2,4,6-trideoxy- β -L-altropyranose + H2O = 5,7-

 $bis (acetylamino) \hbox{--} 3,5,7,9 \hbox{--} tetra de oxy-L-glycero-} \alpha-L-manno-2-nonulo pyranosonic \ acid \ + \ phosphate$

Notes This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1