

## pyruvate kinase

Cat. No. EXWM-3070

Lot. No. (See product label)

## Introduction

**Description** UTP, GTP, CTP, ITP and dATP can also act as donors. Also phosphorylates

hydroxylamine and fluoride in the presence of CO2.

**Synonyms** phosphoenolpyruvate kinase; phosphoenol transphosphorylase

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 2.7.1.40

*CAS No.* 9001-59-6

**Reaction** ATP + pyruvate = ADP + phosphoenolpyruvate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

## Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C $\sim$ -80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com