

butyrate kinase

Cat. No. EXWM-3178

Lot. No. (See product label)

Introduction

Description The enzyme from Clostridium sp. also acts, more slowly, on pentanoate and

propanoate, and on some branched-chain fatty acids (cf. EC 2.7.1.14

sedoheptulokinase).

Product Information

Form Liquid or lyophilized powder

EC Number EC 2.7.2.7

CAS No. 37278-14-1

Reaction ATP + butanoate = ADP + butanoyl phosphate

Notes This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C \sim -80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com