

## acetylsalicylate deacetylase

Cat. No. EXWM-3483

Lot. No. (See product label)

## Introduction

**Description** Not identical with EC 3.1.1.1 (carboxylesterase), EC 3.1.1.2 (arylesterase), EC 3.1.1.7

(acetylcholinesterase) or EC 3.1.1.8 (cholinesterase). The activity of the liver cytosol enzyme is highest with acetyl esters of aryl alcohols, and thioesters are also hydrolysed; the microsomal enzyme also hydrolyses some other negatively charged esters, with highest activity on esters of salicylate with long-

chain alcohols.

**Synonyms** aspirin esterase; aspirin esterase; acetylsalicylic acid esterase; aspirin hydrolase

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 3.1.1.55

*CAS No.* 87348-04-7

**Reaction** acetylsalicylate + H2O = salicylate + acetate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

## Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C $\sim$ -80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1