

## glucuronate-2-sulfatase

Cat. No. EXWM-3751

Lot. No. (See product label)

## Introduction

**Description** Does not act on iduronate 2-sulfate residues (cf. EC 3.1.6.13 iduronate-2-sulfatase)

**Synonyms** glucurono-2-sulfatase

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 3.1.6.18

*CAS No.* 98597-45-6

**Reaction** Hydrolysis of the 2-oulfate groups of the 2-O-sulfo-D-glucuronate residues of

chondroitin sulfate, heparin and heparitin sulfate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

## Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C $\sim$ -80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1