

## geranylgeranyl diphosphate diphosphatase

Cat. No. EXWM-3765

Lot. No. (See product label)

## Introduction

**Synonyms** 

**Description** Involved in the biosynthesis of plaunotol. There are two isoenzymes with different

geranylgeranyl diphosphate phosphatase

ion requirements. Neither require Mg2+ but in addition PII is inhibited by Zn2+, Mn2+ and Co2+. It is not known which isoenzyme is involved in plaunotol

biosynthesis.

•

Product Information

**Form** Liquid or lyophilized powder

**EC Number** EC 3.1.7.5

**Reaction** geranylgeranyl diphosphate + H2O = geranylgeraniol + diphosphate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C∼-80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1