

## endo-α-sialidase

Cat. No. EXWM-3812

Lot. No. (See product label)

## Introduction

**Description** Although the name endo-N-acetylneuraminidase has also been used for this enzyme, this is misleading

since its activity is not restricted to acetylated substrates. An exo-α-sialidase activity is listed as EC

3.2.1.18 exo- $\alpha$ -sialidase. See also EC 4.2.2.15 anhydrosialidase.

**Synonyms** endo-N-acylneuraminidase; endoneuraminidase; endo-N-acetylneuraminidase; poly( $\alpha$ -2,8-sialosyl) endo-N-acylneuraminidase;

acetylneuraminidase; poly( $\alpha$ -2,8-sialoside)  $\alpha$ -2,8-sialosylhydrolase; endosialidase; endo-N

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 3.2.1.129

**CAS No.** 91195-87-8

**Reaction** Endohydrolysis of  $(2\rightarrow 8)-\alpha$ -sialosyl linkages in oligo- or poly(sialic) acids

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

## Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C∼-80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1