

## mucorpepsin

Cat. No. EXWM-4261

Lot. No. (See product label)

## Introduction

**Description** Isolated from the zygomycete fungi Mucor pusillus and M. miehei. The two species variants show 83%

sequence identity and are immunologically crossreactive. In peptidase family A1 (pepsin A family).

Formerly included in EC 3.4.23.6

**Synonyms** Mucor rennin; Mucor aspartic proteinase; Mucor acid proteinase; Mucor acid protease; Mucor miehei

aspartic proteinase; Mucor miehei aspartic protease; Mucor aspartic proteinase; Mucor pusillus emporase;

Fromase 100; Mucor pusillus rennin; Fromase 46TL; Mucor miehei rennin

## **Product Information**

**Form** Liquid or lyophilized powder

**EC Number** EC 3.4.23.23

*CAS No.* 148465-73-0

Reaction Hydrolysis of proteins, favouring hydrophobic residues at P1 and P1'. Clots milk. Does not accept Lys at

P1, and hence does not activate trypsinogen

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

## Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1