

nonpolar-amino-acid-transporting ATPase

Cat. No. EXWM-4659

Lot. No. (See product label)

Introduction

Description Get Nonpolar Amino Acid-Transporting ATPase for advanced amino acid transport research. Creative

Enzymes ensures quality and performance.

Product Information

Form Liquid or lyophilized powder

EC Number EC 7.4.2.2 (Formerly EC 3.6.3.22)

Reaction ATP + H2O + nonpolar amino acidout = ADP + phosphate + nonpolar amino acidin

Notes This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

Storage and Shipping Information

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1