

minus-end-directed kinesin ATPase

Cat. No. EXWM-4707

Lot. No. (See product label)

Introduction

Description Structurally almost identical to EC 3.6.4.3 (microtubule-severing ATPase) but the movement it catalyses is towards the minus end of microtubules.

Product Information

Form Liquid or lyophilized powder

EC Number EC 3.6.4.5

Reaction $\text{ATP} + \text{H}_2\text{O} = \text{ADP} + \text{phosphate}$

Notes This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.