

## minus-end-directed kinesin ATPase

Cat. No. EXWM-4707

Lot. No. (See product label)

### Introduction

#### Description

Structurally almost identical to EC 3.6.4.3 (microtubule-severing ATPase) but the movement it catalyses is towards the minus end of microtubules.

### Product Information

#### Form

Liquid or lyophilized powder

#### EC Number

EC 3.6.4.5

#### Reaction

ATP + H<sub>2</sub>O = ADP + phosphate

#### Notes

This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

### Storage and Shipping Information

#### Storage

Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.