

terpinolene synthase

Cat. No. EXWM-5125

Lot. No. (See product label)

Introduction

Description Requires Mg2+. Mn2+ is less effective and product ratio changes. Forms traces of

other monoterpenoids.

Synonyms ag9; PmeTPS2; LaLIMS_RR

Product Information

Form Liquid or lyophilized powder

EC Number EC 4.2.3.113

Reaction geranyl diphosphate = terpinolene + diphosphate

Notes This item requires custom production and lead time is between 5-9 weeks. We can

custom produce according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1