

β-copaene synthase

Cat. No. EXWM-5140

Lot. No. (See product label)

Introduction

Description Isolated from the fungus Coprinus cinereus. The enzyme also forms (+)- Δ -cadinene, β -cubebene, (+)-

sativene and traces of several other sequiterpenoids. β -Copaene is formed in the presence of Mg2+ but not Mn2+. See EC 4.2.3.13, (+)- Δ -cadinene synthase, EC 4.2.3.128, β -cubebene synthase, and EC

4.2.3.129, (+)-sativene synthase.

Synonyms cop4

Product Information

Form Liquid or lyophilized powder

EC Number EC 4.2.3.127

Reaction (2E,6E)-farnesyl diphosphate = β -copaene + diphosphate

Notes This item requires custom production and lead time is between 5-9 weeks. We can custom produce

according to your specifications.

Storage and Shipping Information

Storage Store it at +4 °C for short term. For long term storage, store it at -20 °C \sim -80 °C.

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1