

## D-alanine-D-serine ligase

Cat. No. EXWM-5752

Lot. No. (See product label)

### Introduction

**Description** The product of this enzyme, D-alanyl-D-serine, can be incorporated into the peptidoglycan pentapeptide instead of the usual D-alanyl-D-alanine dipeptide, which is formed by EC 6.3.2.4, D-alanine-D-alanine ligase. The resulting peptidoglycan does not bind the glycopeptide antibiotics vancomycin and teicoplanin, conferring resistance on the bacteria.

**Synonyms** VanC; VanE; VanG

### Product Information

**Form** Liquid or lyophilized powder

**EC Number** EC 6.3.2.35

**Reaction** D-alanine + D-serine + ATP = D-alanyl-D-serine + ADP + phosphate

**Notes** This item requires custom production and lead time is between 5-9 weeks. We can custom produce according to your specifications.

### Storage and Shipping Information

**Storage** Store it at +4 °C for short term. For long term storage, store it at -20 °C~-80 °C.