

L-Aminoacylase (Crude Enzyme)

Cat. No. NATE-1835

Lot. No. (See product label)

Introduction

Description

This enzyme belongs to the family of hydrolases, those acting on carbon-nitrogen bonds other than peptide bonds, specifically in linear amides. This product with the indicated enzyme activity was briefly purified from engineered E. coli.

Applications

synthesis; medicine;

Synonyms

aminoacylase 1; aminoacylase I; dehydropeptidase II; histozyme; hippuricase; benzamidase; acylase I; hippurase; amido acid deacylase; L-aminoacylase; acylase; aminoacylase; L-amino-acid acylase; α -N-acylaminoacid hydrolase; long acyl amidoacylase; short acyl amidoacylase; ACY1 (gene name); N-acyl-L-amino-acid amidohydrolase

Product Information

Source

E. coli

Appearance

Clear to translucent yellow solution

EC Number

EC 3.5.1.14

CAS No.

9012-37-7

Activity

Undetermined

Reaction

(1) an N-acyl-aliphatic-L-amino acid + H₂O = an aliphatic L-amino acid + a carboxylate (2) an N-acetyl-L-cysteine-S-conjugate + H₂O = an L-cysteine-S-conjugate + acetate

Notes

Since this product needs to be freshly prepared, it will take about 2 weeks after you confirm the order. Each time of the freeze-thawing may cause partial inactivation. Therefore, it should be dispensed as required and stored at -20 °C or lower. With the preservation of the extension of time, the enzyme activity will decline to a certain extent, so the product should be used as soon as possible. This product may have turbidity or precipitation in the production and preservation process, it can be mixed after melting and will not affect the normal use. This product is limited to scientific research use, shall not be used for clinical diagnosis or treatment, shall not be used for food or medicine, shall not be stored in ordinary residential. For your safety and health, please wear an experimental suit and wear disposable gloves.

Usage and Packaging

Package

100ml

Storage and Shipping Information

Storage

at -20 °C or lower, for at least 1 month.