

MoS₂ Nanozymes

Cat. No. NE024

Lot. No. (See product label)

Product Information

Size 10 ml

Parameter 1 mg/ml

**Terence of
Dimension** 0.01

Shape Stripe

**Mimic
Activity** Peroxidase

Description Molybdenum disulfide, as an emerging two-dimensional sheet material like graphene, has been shown to have excellent peroxidase mimetic activity. They can be used to replace biological enzymes in the construction of biosensors.

Application Biosensing, bioimaging, drug delivery, and therapy.

Purity 95%+

**Storage
Condition** Store at 2-8°C in the dark.

Shelf Life 12 months

Usage For Research use only.