

Collagen from Rat Lung

Cat. No. CEFX-006

Lot. No. (See product label)

Introduction

Description Rat lung tissue was taken to prepare collagen by pepsin-acid extraction, salt

precipitation, re-solubilization in acetic acid, and dialysis and re-precipitation in

neutral phosphate buffer. Collagen was re-solubilized and lyophilized.

Electrophoresis in SDS polyacrylamide showed that the collagen was of high purity

with typical migration of $\gamma,\,\beta(2),\,\alpha 1,$ and $\alpha 2$ bands.

Product Information

Species Rat

Source Rat Lung

Form Powder

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1