

## Natural Liquid Bandage

Cat. No. BPEL-016

Lot. No. (See product label)

### Introduction

#### Description

This product is composed of complex phage, alginate, proline, sodium hyaluronate and so on. Complex phage can be quickly sterilized, and will not cause bacterial imbalance. Alginate forms a gel layer on the wound surface to resist bacterial erosion and keep the wound moist, absorbing the secreted pus and accelerating healing. Proline helps the body break down proteins that are used to create healthy cells in the body. Important for maintaining healthy skin and connective tissue, especially in traumatized areas. Sodium hyaluronate supports the proliferation and migration of epidermal fibroblasts to help maintain normal skin hydration.

### Product Information

#### Appearance

Liquid

#### Function

Antibacterial and anti-inflammatory, pain relief, promote wound healing. Used for surgical or traumatic (e.g. broken tail, standing ear, cesarean section, lacerations, etc.) wound treatment. For the treatment of various skin wound infections as well as the treatment of burns.

### Storage and Shipping Information

#### Storage

Store in a cool, dry place. Shelf life: 12 months.