

## Phosphatase Inhibitor Cocktail B

Cat. No. CEI-0786

Lot. No. (See product label)

### Introduction

**Description** Reconstitute in 1ml water to make a solution with 200mM imidazole, 100mM sodium fluoride, 115mM sodium molybdate, 100mM sodium orthovanadate and 400mM sodium tartrate. For the inhibition of protein tyrosine phosphatases, acid- and alkaline-phosphatases.

### Product Information

**Appearance** Lyophilized solid

**Targets** Phosphatase