

## Phosphatase Inhibitor Cocktail II

Cat. No. CEI-0790

Lot. No. (See product label)

### Introduction

**Description** 100X concentrate. This cocktail will inhibit broad spectrum of serine/threonine and tyrosine phosphatases and is suitable for cell and tissue lysis. Contains: 100mM sodium fluoride, 50mM sodium orthovanadate, 50mM sodium pyrophosphate, and 100mM β-glycerophosphate.

### Product Information

**Appearance** Liquid

**Targets** Phosphatase

### Storage and Shipping Information

**Storage** Store at -20°C.

**Shipping** Not regulated for shipping

**Conditions**