

Phosphatase Inhibitor Cocktail II

Cat. No. CEI-0790

Lot. No. (See product label)

Introduction

Description

100X concentrate. This cocktail will inhibit broad spectrum of serine/threonine and tyrosine phosphatases and is suitable for cell and tissue lysis. Contains: 100mM sodium fluoride, 50mM sodium orthovanadate, 50mM sodium pyrophosphate, and 100mM b-glycerophosphate.

Product Information

Appearance

Liquid

Targets

Phosphatase

Storage and Shipping Information

Storage

Store at -20°C.

Shipping Conditions

Not regulated for shipping