

Protease Inhibitor Cocktail I

Cat. No. CEI-0794

Lot. No. (See product label)

Introduction

Description Forms 100x solution when reconstituted in 1ml water. A 1x solution contains

500uM AEBSF.HCl, 150nM aprotinin, 1uM E-64, 0.5mM EDTA disodium salt, and 1uM

leupeptin hemisulfate.

Product Information

Appearance Lyophilized solid

Targets Protease

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1