

Protease Inhibitor Cocktail III

Cat. No. CEI-0797

Lot. No. (See product label)

Introduction

Description 1ml in DMSO. One unit contains 100mM AEBSF.HCl, 80uM aprotinin, 5mM bestatin, 1.5mM E-64, 2mM

leupeptin hemisulfate, and 1mM pepstatin A. For use with mammalian cell and tissue extracts.

Product Information

Appearance Liquid

Targets Protease

Storage and Shipping Information

Shipping

Not regulated for shipping

Conditions

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1