

Protease Inhibitor Cocktail V, EDTA-Free

Cat. No. CEI-0801

Lot. No. (See product label)

Introduction

Description For inhibition of cysteine and serine proteases, but not not metalloproteases.

Reconstitute in 1ml of water for a 100X solution. A 1X solution contains $500\mu M$

AEBSF.HCl, 150nM aprotinin, 1μ M E-64 and 1μ M leupeptin hemisulfate.

Product Information

Appearance Lyophilized solid

Targets Protease

Tel: 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1