

## Liberase MNP-S GMP Grade

Cat. No. NATE-0996

Lot. No. (See product label)

### Introduction

**Applications** Use Liberase MNP-S, GMP Grade for in vitro enzymatic dissociation of tissue, to isolate single cells from a broad range of tissue types. In particular intended for use in isolation of tissue based stem cells and chondrocytes from cartilage.

**Synonyms** Liberase; Liberase MNP-S

### Product Information

**Appearance** White lyophilized cake

**Activity** 19,140-37,000U/bottle

### Storage and Shipping Information

**Stability** At -15 to -25°C within specification range for 24 months.