

## Lysozyme 25A from *Streptococcus pneumoniae*, Recombinant

Cat. No. NATE-1456

Lot. No. (See product label)

### Introduction

#### Description

Lysozymes, also known as muramidase or N-acetylmuramide glycanhydrolase, are glycoside hydrolases. These are enzymes (EC 3.2.1.17) that damage bacterial cell walls by catalyzing hydrolysis of 1,4-beta-linkages between N-acetylmuramic acid and N-acetyl-D-glucosamine residues in a peptidoglycan and between N-acetyl-D-glucosamine residues in chitodextrins. Lysozyme is abundant in a number of secretions, such as tears, saliva, human milk, and mucus. It is also present in cytoplasmic granules of the macrophages and the polymorphonuclear neutrophils (PMNs). Large amounts of lysozyme can be found in egg white. C-type lysozymes are closely related to alpha-lactalbumin in sequence and structure, making them part of the same family. In humans, the lysozyme enzyme is encoded by the LYZ gene.

#### Synonyms

muramidase; globulin G; mucopeptide glucohydrolase; globulin G1; N,O-diacetylmuramidase; lysozyme g; L-7001; 1,4-N-acetylmuramidase; mucopeptide N-acetylmuramoylhydrolase; PR1-lysozyme; lysozyme; LYZ; LZM; EC 3.2.1.17; 9001-63-2

### Product Information

<b>Species</b>	<i>Streptococcus pneumoniae</i>
<b>Source</b>	<i>E. coli</i>
<b>Form</b>	35 mM NaHepes buffer, pH 7.5, 750 mM NaCl, 200 mM imidazol, 3.5 mM CaCl <sub>2</sub> , 0.02% sodium azide and 25% (v/v) glycerol
<b>EC Number</b>	EC 3.2.1.17
<b>CAS No.</b>	9001-63-2
<b>Molecular Weight</b>	26.3 kDa
<b>Purity</b>	>90% by SDS-PAGE
<b>Concentration</b>	1 mg/mL
<b>Optimum pH</b>	5
<b>Optimum temperature</b>	37 °C
<b>Specificity</b>	4-nitrophenyl-(GlcNAc) <sub>5</sub>

### Storage and Shipping Information

<b>Storage</b>	This enzyme is shipped at room temperature but should be stored at -20 °C.
----------------	--