

## Chymotrypsin-polyethylene glycol

*Cat. No.* NATE-1581

*Lot. No.* (See product label)

### **Introduction**

**Description** Enzyme PEG conjugate was made from by the covalent attachment of PEG 5000 to the enzyme molecule via its amino groups. Conjugate is provided in lyophilized powder form.

**Synonyms** Chymotrypsin PEG Conjugate; Chymotrypsin-PEG; PEG-Chymotrypsin; Chymotrypsin-polyethylene glycol

### **Product Information**

**Form** Lyophilized powder

### **Storage and Shipping Information**

**Storage** Store at -20°C