

# 3,3'-Diaminobenzidine tetrahydrochloride

Cat. No. CSUB-0419

Lot. No. (See product label)

### Introduction

ApplicationsReagent for spectrophotometric determination of selenium.

**Synonyms** 3, 3'; 4, 4'-Biphenyltetramine tetrahydrochloride; 3, 3'; 4, 4'-Tetraaminobiphenyl

tetrahydrochloride; DAB

#### **Product Information**

**Form** tablet

**EC Number** 231-018-9

*CAS No.* 7411-49-6

Molecular Weight 360.11

**Specificity** Peroxidase substrate suitable for use in immunoblotting and immunohistological

staining procedures. This substrate produces an insoluble end product that is

brown in color and can be observed visually.

**Substrates** Peroxidase Substrates

## **Usage and Packaging**

**Reconstitution** Dissolve 1 tablet in 15 ml of Tris-buffered saline, pH 7.6. Add 12 µl of fresh 30%

hydrogen peroxide. DAB and hydrogen peroxide concentration may be adjusted to

1/1

suit individual applications.

## Storage and Shipping Information

*Storage* -20°C

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com