

## 3,3',5,5'-Tetramethylbenzidine Liquid Substrate, Super Slow, for ELISA

Cat. No. CSUB-0461

Lot. No. (See product label)

Introduction

**Applications** Contains TMB in a mildly acidic buffer. The substrate is supplied as a one

component ready to use solution. Rate kinetics are approx. 60% slower than

traditional TMB formulations. Not recommended for membrane or

immunohistochemical applications that require a precipitating reaction product.

**Product Information** 

**Substrates** Peroxidase Substrates

Storage and Shipping Information

**Storage** 2-8°C

**Shipping Conditions** wet ice

**Tel:** 1-631-562-8517 1-516-512-3133 **Email:** info@creative-enzymes.com

1/1