

## Urea solution

Cat. No. CSUB-0608

Lot. No. (See product label)

### Introduction

**Applications** Recommended for the detection of urease activity.

**Synonyms** Carbamide; Carbonyldiamide

### Product Information

**CAS No.** 57-13-6

**Molecular Weight** 60.06

**Concentration** 40 % (w/v) in H<sub>2</sub>O

**Substrates** Urease

### Usage and Packaging

**Preparation Instructions** Warm the refrigerated 40% urea solution to room temperature. Aseptically add 5 mL to 95 mL of sterile molten Urea Agar Base (Christensen). Mix well and dispense in sterile tubes.

### Storage and Shipping Information

**Storage** 2-8°C